

Product Datasheet

NT-proBNP Antigen, Recombinant >95% LA313

Product Name	NT-proBNP Antigen, Recombinant >95%
Catalog Number	LA313
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0 ± 0.2
Matched Pairs	HM071
Expression System	E. coli
Preservatives	0.09% Sodium Azide