

Product Datasheet

Treponema pallidum (TpN15) Antigen, Recombinant >95% LA314

Product Name	Treponema pallidum (TpN15) Antigen, Recombinant >95%
Catalog Number	LA314
Tested Applications	ELISA, WB
Alternative Names	Syphilis
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	>95%
Buffer	25 mM Tris-HCl, pH 8.3 with 0.1% SDS
Specificity	TpN15
Molecular Weight	40.7 kDa
Expression System	E. coli