

## Product Datasheet

### Human CD141 ELISA Set 851.690

<b>Product Name</b>	Human CD141 ELISA Set
<b>Catalog Number</b>	851.690
<b>Description</b>	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
<b>Range</b>	0.62-20 ng/ml
<b>Sensitivity</b>	0.31 ng/ml
<b>Alternative Names</b>	Thrombomodulin
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Specificity</b>	Recognizes both natural and recombinant human sCD141 (Thrombomodulin)
<b>Sample Types</b>	Serum, Plasma, Cell culture supernatant