

Product Datasheet

Norovirus genogroup II Antigen, Recombinant >95% LA321

Product Name	Norovirus genogroup II Antigen, Recombinant >95%
Catalog Number	LA321
Tested Applications	ELISA, Control
Immunogen	Norovirus GII.4 VLP
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (Full length VP1 peptide with His-tag at N-terminus)
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline with 50 mM Argine
Specificity	Norovirus GII
Matched Pairs	HM638 or HM885
Expression System	E. coli
Preservatives	0.05% Sodium Azide