

Product Datasheet

Human CD141 ELISA Set 851.690.005

Product Name	Human CD141 ELISA Set
Catalog Number	851.690.005
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	0.62-20 ng/ml
Sensitivity	0.31 ng/ml
Tested Applications	ELISA
Alternative Names	Thrombomodulin
Brand	Diaclone
Storage	Store at +4°C
Specificity	Recognizes both natural and recombinant human sCD141 (Thrombomodulin)
Sample Types	Serum, Plasma, Cell culture supernatant