

Product Datasheet

Herpes Simplex Virus 2 (HSV 2-gG) Antigen, Recombinant >90% LA328

Product Name	Herpes Simplex Virus 2 (HSV 2-gG) Antigen, Recombinant >90%
Catalog Number	LA328
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	≤-70°C
Purity	>90%
Source	Pichia pastoris
Buffer	50 mM Tris, 5mM DTT, 1mM EDTA, 0.1% Triton X-100, pH 8.0
Specificity	gG
Expression System	Pichia pastoris