

Product Datasheet

Troponin I, Cardiac Antigen, Recombinant >95% LA369L

Product Name	Troponin I, Cardiac Antigen, Recombinant >95%
Catalog Number	LA369L
Tested Applications	WB, ELISA
Brand	EastCoast Bio
Form/Appearance	Lyophilized Purified Recombinant Antigen (His-tag)
Storage	Short Term (< 1 week): 2-8°C. Long Term: ≤-20°C. After reconstitution: -20°C. Avoid repeated freezing and thawing.
Note	Reconstitute with water to 4 mg/ml.
Purity	>95% by SDS-PAGE
Molecular Weight	26 kDa
Expression System	E. coli