

## Product Datasheet

### Creatine Kinase BB (CK-BB) Antigen, Recombinant >95% LA374

<b>Product Name</b>	Creatine Kinase BB (CK-BB) Antigen, Recombinant >95%
<b>Catalog Number</b>	LA374
<b>Tested Applications</b>	WB, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen (His-tag)
<b>Storage</b>	Short Term ( $\leq 1$ week): 2-8°C. Long Term: $\leq -20^{\circ}\text{C}$ . Avoid repeated freezing and thawing.
<b>Note</b>	Amino Acid Sequence: MPFSNSHNALKLRFPAEDEFPDLSAHNNHMAKVLTPELYAELRAKSTPSGFTLDDVIQTGVDNPGHPYIMTVGCVAGDE
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Molecular Weight</b>	44 kDa
<b>Expression System</b>	E. coli