

# Product Datasheet

## C-Reactive Protein (CRP) Antigen, Native >95% LA378

Product Name	C-Reactive Protein (CRP) Antigen, Native >95%
Catalog Number	LA378
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Liquid, 0.2 µm filtered
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Tested negative for HBsAg, HBV, HCV, and HIV 1 & 2 by FDA approved methods. However material should be handled as potentially infectious.
Purity	>95%
Source	Human Fluids
Buffer	0.1M Tris, pH 7.5 with 0.2M NaCl and 2 mM CaCl <sub>2</sub>
Preservatives	0.1% Sodium Azide