

Product Datasheet

Ebola Virus Antigen, Recombinant >95% LA382

Product Name	Ebola Virus Antigen, Recombinant >95%
Catalog Number	LA382
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Arginin
Specificity	Ebola Virus (Zaire) IgG and IgM
Expression System	E. coli
Preservatives	0.02% Sodium Azide