

Product Datasheet

UR-144-BSA Antigen LA385

Product Name	UR-144-BSA Antigen
Catalog Number	LA385
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA Conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	20 mM Phosphate Buffered Saline, pH 7.5
Matched Pairs	HM595 or HM596 (detection)
Preservatives	0.01% Sodium Azide