

# Product Datasheet

## Human IL-17A ELISA Set 851.720

Product Name	Human IL-17A ELISA Set
Catalog Number	851.720
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	3.125-100 pg/ml
Sensitivity	3 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-17A
Sample Types	Serum, Plasma, Cell culture supernatant