

# Product Datasheet

## Human IL-17A ELISA Set 851.720.005

|                     |  |
|---------------------|--|
| Product Name        | Human IL-17A ELISA Set   |
| Catalog Number      | 851.720.005  |
| Description         | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range               | 3.125-100 pg/ml  |
| Sensitivity         | 3 pg/ml  |
| Tested Applications | ELISA  |
| Alternative Names   | NULL   |
| Brand               | Diaclone   |
| Storage             | Store at +4°C  |
| Cross Reactivity    | No cross reactivity with other human cytokines   |
| Specificity         | Recognizes both natural and recombinant human IL-17A   |
| Sample Types        | Serum, Plasma, Cell culture supernatant  |