

Product Datasheet

Helicobacter pylori (OMP) Antigen, Recombinant LA396

Product Name	Helicobacter pylori (OMP) Antigen, Recombinant
Catalog Number	LA396
Species/Host	E. coli
Tested Applications	LF, ELISA, TIA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	2-8°C
Purity	≥95%
Buffer	Phosphate buffer with additives and salt, pH 7.0 ± 0.2
Matched Pairs	HM629
Expression System	E. coli
Preservatives	0.09% Sodium Azide