

## Product Datasheet

### Rotavirus (VP6) Antigen, Recombinant >90% LA397

**Product Name** Rotavirus (VP6) Antigen, Recombinant >90%

**Catalog Number** LA397

**Tested Applications** ELISA, LF

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen

**Storage** Short Term: 2-8°C. Long Term: -20°C

**Purity** >90%

**Buffer** Tris Buffer with additives and salts, pH 7.5

**Matched Pairs** HM630

**Expression System** E. coli

**Preservatives** 0.09% Sodium Azide