

Product Datasheet

Calprotectin Antigen, Recombinant >95% LA399

Product Name Calprotectin Antigen, Recombinant >95%

Catalog Number LA399

Tested Applications LF, ELISA, TIA

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage 2-8°C

Purity >95%

Source E.coli

Buffer Tris Buffer, pH 8.0

Molecular Weight 24 kDa

Matched Pairs HM632

Expression System E.coli

Preservatives 0.09% Sodium Azide