

## Product Datasheet

### Calprotectin Antigen, Recombinant >95% LA399

**Product Name** Calprotectin Antigen, Recombinant >95%

**Catalog Number** LA399

**Tested Applications** LF, ELISA, TIA

**Brand** EastCoast Bio

**Form/Appearance** Purified Recombinant Antigen (His-tag)

**Storage** 2-8°C

**Purity** >95%

**Source** E.coli

**Buffer** Tris Buffer, pH 8.0

**Molecular Weight** 24 kDa

**Matched Pairs** HM632

**Expression System** E.coli

**Preservatives** 0.09% Sodium Azide