

# Product Datasheet

## Zika Virus Antigen, Native LA418

Product Name	Zika Virus Antigen, Native
Catalog Number	LA418
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Native Antigen
Storage	≤-70°C
Note	Although material has been inactivated and tested negative with infectivity challenge, as with all biological reagents, it should be treated as potentially infectious.
Source	E6 Cells/Strain MR 766
Buffer	50 mM Glycine with 140 mM NaCl, pH 9.0