

Product Datasheet

Escherichia coli (E. coli) O157 Antigen, Native LA420

Product Name	Escherichia coli (E. coli) O157 Antigen, Native
Catalog Number	LA420
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Native Inactivated Antigen
Storage	Short Term (< 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM649
Preservatives	0.09% Sodium Azide