

Product Datasheet

Dengue Virus (DENV-Env Serotype 2) Antigen, Recombinant >90% LA421

Product Name	Dengue Virus (DENV-Env Serotype 2) Antigen, Recombinant >90%
Catalog Number	LA421
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>90%
Buffer	20 mM Phosphate Buffered Saline
Expression System	E. coli
Preservatives	0.01% ProClin 300