

Product Datasheet

Salmonella typhi Antigen, Native LA424

Product Name	Salmonella typhi Antigen, Native
Catalog Number	LA424
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Inactivated Native Antigen
Storage	Short Term (< 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Source	S. typhi culture
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM174, HM037, HM038, HM039, or HM040
Preservatives	0.09% Sodium Azide