

## Product Datasheet

### Salmonella paratyphi A Antigen, Native LA425

<b>Product Name</b>	Salmonella paratyphi A Antigen, Native
<b>Catalog Number</b>	LA425
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Inactivated Native Antigen
<b>Storage</b>	Short Term (< 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Source</b>	S. paratyphi A culture
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Matched Pairs</b>	HM175, HM514, HM037, HM038, HM039, or HM040
<b>Preservatives</b>	0.09% Sodium Azide