

Product Datasheet

Salmonella paratyphi A Antigen, Native LA425

Product Name	Salmonella paratyphi A Antigen, Native
Catalog Number	LA425
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Inactivated Native Antigen
Storage	Short Term (< 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Source	S. paratyphi A culture
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 ± 0.2
Matched Pairs	HM175, HM514, HM037, HM038, HM039, or HM040
Preservatives	0.09% Sodium Azide