

Product Datasheet

Aldosterone-BSA Antigen LA427

Product Name	Aldosterone-BSA Antigen
Catalog Number	LA427
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM651, HM652, HM653 or HM783 (detection)
Preservatives	0.01% Sodium Azide