

Product Datasheet

Estradiol (E2)-6-BSA Antigen LA429

Product Name	Estradiol (E2)-6-BSA Antigen
Catalog Number	LA429
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: ≤-20°C. Avoid repeated freezing and thawing.
Purity	≥95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 ± 0.2
Matched Pairs	HM654 (detection)
Preservatives	0.01% Sodium Azide