

Product Datasheet

Corticosterone-3-BSA Antigen LA431

Product Name	Corticosterone-3-BSA Antigen
Catalog Number	LA431
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: ≤-20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM655 or HM656 (detection)
Preservatives	0.05% Sodium Azide