

Product Datasheet

Human CD178 ELISA Set 851.730.005

Product Name	Human CD178 ELISA Set
Catalog Number	851.730.005
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	62.5-2000 pg/ml
Sensitivity	12 pg/ml
Tested Applications	ELISA
Alternative Names	FasL, CD95L
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes both natural and recombinant soluble human Fas-Ligand
Sample Types	Serum, Plasma, Cell culture supernatant