

Product Datasheet

Hepatitis B Core Antigen (HBcAg) Antigen, Recombinant >90% LA439

Product Name	Hepatitis B Core Antigen (HBcAg) Antigen, Recombinant >90%
Catalog Number	LA439
Tested Applications	WB, ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
Purity	>90%
Buffer	20 mM Tris-HCl, pH 8.0 with 1 mM EDTA
Molecular Weight	15 kDa
Expression System	E. coli
Preservatives	0.1% Proclin 300