

Product Datasheet

Hepatitis Be Antigen (HBeAg), Recombinant >90% LA440

Product Name	Hepatitis Be Antigen (HBeAg), Recombinant >90%
Catalog Number	LA440
Tested Applications	WB, ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen (His-tag)
Storage	Short Term: 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
Purity	>90%
Buffer	20 mM Tris-HCl, pH 8.0 with 1 mM EDTA
Molecular Weight	14 kDa
Expression System	E. coli
Preservatives	0.1% Proclin 300