

## Product Datasheet

### Carcinoembryonic Antigen (CEA), Recombinant >95% LA444

<b>Product Name</b>	Carcinoembryonic Antigen (CEA), Recombinant >95%
<b>Catalog Number</b>	LA444
<b>Tested Applications</b>	Control
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	≤-20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	>95% by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Molecular Weight</b>	100-180 kDa
<b>Expression System</b>	Mammalian Cell
<b>Preservatives</b>	0.02% Proclin 300