

Product Datasheet

Lipoprotein a (Lp-a) Antigen, Native >90% LA447

Product Name	Lipoprotein a (Lp-a) Antigen, Native >90%
Catalog Number	LA447
Description	Purified Native Antigen
Tested Applications	Immunoassays, Control
Brand	EastCoast Bio
Form/Appearance	Liquid
Storage	≤-20°C. Avoid repeated freezing and thawing.
Purity	≥90% by SDS-PAGE
Source	Human Serum
Buffer	0.02M PBS, pH 7.4
Preservatives	0.05% ProClin 300