

Product Datasheet

Inhibin B (INHB) Antigen, Recombinant >95% LA448

Product Name Inhibin B (INHB) Antigen, Recombinant >95%

Catalog Number LA448

Tested Applications ELISA, CLIA, LF

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage ≤ -20° C

Purity >95% by SDS-PAGE

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Molecular Weight 46 kDa

Expression System Mammalian Cell

Preservatives 0.02% Proclin 300