

Product Datasheet

Gastrin 17 (G-17) Antigen, Synthetic >95% LA451

Product Name Gastrin 17 (G-17) Antigen, Synthetic >95%

Catalog Number LA451

Tested Applications Control

Brand EastCoast Bio

Form/Appearance Purified Synthetic Antigen

Storage ≤ -20°C. Avoid repeated freezing and thawing.

Purity >95% by SDS-PAGE

Source Synthetic

Buffer DMSO

Matched Pairs HM774, HM775, HM776