

Product Datasheet

Helicobacter pylori (CagA) Antigen, Recombinant $\geq 90\%$ LA460

Product Name	Helicobacter pylori (CagA) Antigen, Recombinant $\geq 90\%$
Catalog Number	LA460
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Chimeric Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	$\geq 90\%$
Buffer	10 mM Phosphate Buffer, pH 7.4 \pm 0.2 with 5mM EDTA, 0.01% SDS
Specificity	H. pylori CagA & GroEL
Molecular Weight	98 kDa
Matched Pairs	LA460 (self-pairs)
Expression System	E. coli