

Product Datasheet

Buprenorphine (BUP)-BSA Antigen LA463

Product Name	Buprenorphine (BUP)-BSA Antigen
Catalog Number	LA463
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	Phosphate Buffered Saline, pH 7.2
Preservatives	0.02% Sodium Azide