

Product Datasheet

Benzoylecgonine (BE)-BTG Antigen LA466

Product Name	Benzoylecgonine (BE)-BTG Antigen
Catalog Number	LA466
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BTG conjugate
Storage	Short Term (≤ 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	$\geq 90\%$ by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM777 (detection)
Preservatives	0.01% Sodium Azide