

Product Datasheet

Astrovirus Antigen, Recombinant LA470

Product Name	Astrovirus Antigen, Recombinant
Catalog Number	LA470
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	TBS, pH 8.0
Matched Pairs	HM803, HM804
Expression System	E. coli
Preservatives	0.09% Sodium Azide