

Product Datasheet

Enterovirus (VP1) Antigen, Recombinant LA471

Product Name	Enterovirus (VP1) Antigen, Recombinant
Catalog Number	LA471
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Recombinant Antigen
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	TBS, pH 7.4
Matched Pairs	JP134
Expression System	E. coli
Preservatives	0.09% Sodium Azide