

Product Datasheet

Campylobacter Antigen, Native LA473

Product Name	Campylobacter Antigen, Native
Catalog Number	LA473
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Native Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Source	Campylobacter culture
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM805, HM806
Preservatives	0.09% Sodium Azide