

## Product Datasheet

### Campylobacter jejuni Antigen, Native LA474

<b>Product Name</b>	Campylobacter jejuni Antigen, Native
<b>Catalog Number</b>	LA474
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Native Antigen
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Source</b>	C. jejuni culture
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Matched Pairs</b>	HM805, HM806
<b>Preservatives</b>	0.09% Sodium Azide