

## Product Datasheet

### Lactoferrin Antigen, Native LA478

<b>Product Name</b>	Lactoferrin Antigen, Native
<b>Catalog Number</b>	LA478
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Native Antigen
<b>Storage</b>	2-8°C
<b>Purity</b>	≥ 95 %
<b>Source</b>	Human Milk
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.4
<b>Matched Pairs</b>	HM810, HM811
<b>Preservatives</b>	0.09% Sodium Azide