

Product Datasheet

Human IL-17F ELISA Set 851.740.015

Product Name	Human IL-17F ELISA Set
Catalog Number	851.740.015
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	15.6-500 pg/ml
Sensitivity	6.0 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-17F
Sample Types	Serum, Plasma, Cell culture supernatant