

## Product Datasheet

### Elastase CELA3B Antigen, Recombinant LA481

<b>Product Name</b>	Elastase CELA3B Antigen, Recombinant
<b>Catalog Number</b>	LA481
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Recombinant Antigen
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	Tris buffer with additives and salts, pH 8.0
<b>Matched Pairs</b>	JP137, JP138
<b>Expression System</b>	Mammalian Cell
<b>Preservatives</b>	0.09% Sodium Azide