

Product Datasheet

Ebola Virus (Zaire-NP) Antigen, Recombinant >95% LA489

Product Name	Ebola Virus (Zaire-NP) Antigen, Recombinant >95%
Catalog Number	LA489
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	-20°C. Avoid repeated freezing and thawing.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Ebola Virus Zaire nucleoprotein
Molecular Weight	15 kDa
Expression System	E.coli
Preservatives	0.02% Sodium Azide