

Product Datasheet

Serum Amyloid A (SAA), Feline Antigen, Native LA493

Product Name	Serum Amyloid A (SAA), Feline Antigen, Native
Catalog Number	LA493
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Native Antigen
Storage	Short Term: -20°C. Long Term: -70°C
Purity	>90%
Source	Feline Serum
Buffer	6M Urea, 100 mM Tris-HCl, 1 mM EDTA, pH 7.5
Molecular Weight	12 kDa
Matched Pairs	HM702, HM703