

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >90% LA497

| | |
|----------------------------|---|
| Product Name | Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >90% |
| Catalog Number | LA497 |
| Tested Applications | ELISA, LF |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Recombinant Antigen (His-tag) |
| Storage | 2-8°C. Do not freeze. |
| Purity | >90% |
| Buffer | 20 mM Hepes, 300 mM NaCl, 102 mM Imidazole, 0.2% SDS |
| Specificity | gp41 |
| Molecular Weight | 18 kDa |
| Expression System | E. coli |
| Preservatives | 0.05% Proclin 300 |