

# Product Datasheet

## Human IL-17A/F ELISA Set 851.750.005

Product Name	Human IL-17A/F ELISA Set
Catalog Number	851.750.005
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	15.6-500 pg/ml
Sensitivity	6.6 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-17A/F
Sample Types	Serum, Plasma, Cell culture supernatant