

Product Datasheet

Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >90% LA498

Product Name Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >90%

Catalog Number LA498

Tested Applications ELISA, LF

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen (His-tag)

Storage 2-8°C. Do not freeze.

Purity >90%

Buffer 20 mM HEPES, 300 mM NaCl, 0.2% SDS

Specificity gp36

Molecular Weight 20 kDa

Expression System E. coli

Preservatives 0.05% Proclin 300