

## Product Datasheet

### Trypanosoma cruzi Chimeric Antigen,, Recombinant >90% LA499L

<b>Product Name</b>	Trypanosoma cruzi Chimeric Antigen,, Recombinant >90%
<b>Catalog Number</b>	LA499L
<b>Tested Applications</b>	ELISA, LF
<b>Alternative Names</b>	Chagas disease
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Lyophilized Recombinant Antigen (His-tag)
<b>Storage</b>	2-8°C. ≤ -20°C after reconstitution. Avoid repeated freezing and thawing.
<b>Note</b>	Reconstitute with 500 µl of water and gentle shaking
<b>Purity</b>	>90%
<b>Buffer</b>	20 mM Tris, 100 mM NaCl
<b>Molecular Weight</b>	9 kDa
<b>Expression System</b>	E. coli
<b>Preservatives</b>	0.05% Proclin 300