

# Product Datasheet

## Hepatitis B Surface Antigen (HBsAg-subtype Ay), Native >99% LA501

Product Name	Hepatitis B Surface Antigen (HBsAg-subtype Ay), Native >99%
Catalog Number	LA501
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified & 0.2 µm filtered Native Antigen
Storage	2-8°C
Note	Tested negative for HIV 1/2 and HCV, however as with all biological material it should be handled as potentially infectious.
Purity	>99%
Source	Human Plasma
Buffer	25 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl, 5% sucrose
Preservatives	0.1% Sodium Azide